



Analysis of 9 days dataset – 4th to 12th May

- Expect channels 0 to 22 to be similar (new PMTs), tested by overlapping channels 6 and 7
- Repeated the above for old and new PMTs connected to calo 18 – ch 8 (new) and ch 28 (old)
- A study of channel 8 (new PMT) and 28 (old PMT) connected to calo 18.
- Displayed A1 and A2 on the same plot and A2:A1 separately
- Just mentioned max. range of variation of A1 and A2 (straight line fit did not make sense especially for channel 8).
- Fitted A2:A1 with a straight line as before.
- Finally found an abrupt change in temperature data, so split this dataset into two parts of about 4.5 days each




Temperature Plots - A1 with optical table temperature



Something changed abruptly in the laser hut temperature, so used optical table temperature for analysis. Also divided the dataset in 2 sets – 0 to 5700s and 6000s to 12000s.






Channels 6 and 7 - comparison





Purpose – expected that all new PMTs show similar behaviour (these are comparable). Checked channel 6 and 7. Normalized with the mean of first 9 bins.

Calo 18 – Channels **28** and **8** - comparison




Should show different behaviour (as the PMTs are different and old PMTs use a PMMA return fibre and Silica forward fiber whereas the new one uses only Silica). This is evident here.



Temperature studies of channel 8 – Linear fit




New PMT – Silica fiber both
forward and return directions







Temperature studies of channel 28 – Linear fit





Old PMT – PMMA fiber in the return direction. Thus
A2 = Silica forward + PMMA return




Investigation for the two sets (or parts) separately






Calo 18 – Channels **28** and **8** – comparison – Set I



Should show different behaviour (as the PMTs are different and old PMTs use a PMMA return fibre and Silica forward fiber whereas the new one uses only Silica). This is evident here.





Calo 18 – Channels **28** and **8** – comparison – Set II




Should show different behaviour (as the PMTs are different and old PMTs use a PMMA return fibre and Silica forward fiber whereas the new one uses only Silica). This is evident here.



Temperature studies of channel 8 – Set 1




New PMT – Silica fiber both
forward and return directions



Temperature studies of channel 28 – Set I



Old PMT – PMMA fiber in the return direction. Thus A2 = Silica forward + PMMA return







Conclusions



- Variation of A1 with temperature shows a negative correlation of PMT gain with temperature (as expected)
- A2:A1 is almost constant with time for the (for longer time ranges some the variation appears especially for the old PMTs)
- Channels 0 to 23 i.e. the new PMTs (refer slide 9) show a comparatively stable A2:A1 drop with temperature (it ranges from $\sim 0.1\%$ to 0.45%) which is much better than the old PMTs
- Old PMTs show a much larger change with temperature (probably expected) – besides the return fiber of all PMTs is a PMMA, which has shown to display a greater loss.

Back up slides


Calo 18 - A1 and A2 overlaid and A2:A1





Drop in gain for rest channels – 0 to 23




Results obtained by
linear fits of each
channel.





Temperature studies of channel 8




Used quadratic fit for A1 and A2
– gives better chi2. A cubic fir for
A2:A1






Temperature studies of channel 28



Used higher order fits –
gives better chi2 than linear
fit. A cubic fir for A2 and
A2:A1



Comparison of Old and New PMTs – Gain drop with temp



Results obtained by linear fits of each calorimeter.