

WFD timing

Nandita Raha – INFN Pisa


21 – Mar – 2019

How to get time stability within fill with picosecond precision

1. Check Δt between consecutive pulses
 - consecutive pulses belong to different fills \rightarrow subtract SYNC
 - shift between consecutive pulses is set in $2.5 \mu\text{s}$ steps \rightarrow must collect large statistics
2. Use Short Time Double Pulse
 - STDP (0-100 ns) is normally taken in 4 or 5 moments during fill at distances $\sim 100\text{-}150 \mu\text{s}$ \rightarrow this is done to increase statistics
 - comparison of the Δt set by the Delay Generator with the one measured by the WFD

Δt between consecutive pulses


- The time between consecutive pulses is set in 2.5 μs steps \rightarrow many events required to reach ps precision – used 60 hr dataset
- Preliminary result: **Bin size: Y – axis 12.5 ns, X- axis 5 μs**



Δt between consecutive pulses


- Combined the Δt between consecutive pulses (both first pulse and second pulse in a fill) in a single plot

Bin size: Y – axis 12.5 ps, X- axis 5 μs




Δt between consecutive pulses

Δt consecutive pulses (ns) for 1 bin – three consecutive bins shown



Δt between consecutive pulses


Profiles with linear fit – Mean 2499.891 ± 0.00045 ns




Δt between consecutive pulses

Combined the Δt between consecutive pulses (both first pulse and second pulse in a fill) in a single plot

Profiles with linear fit (pol1) for entire range– Mean 2499891.35 ± 1.26 ps




Δt between consecutive pulses – all 1296 xtal




Obtained by Gaussian fit of Y projection with $\sigma/\sqrt{\text{NDF}}$ as errors.
Error bars of ~ 10 ps

Δt between consecutive pulses – Average for all Lasers



This is the average Δt of all crystals of a single laser (from previous slide).


Short Time Double Pulse




- measured Δt stable at <100 ps level
- shift due to different light path of the two lasers
- residual trend at 5 ps level

Short Time Double Pulse

Difference from Δt_{set}




- measured Δt stable at $<100 \text{ ps}$ level
- shift due to different light path of the two lasers
- residual trend at 5 ps level




- $\Delta(\Delta t)$ stable within 10 ps in 5 different time slices within fill
- This technique is more precise than InFill pulses, but less points within fill available

$(\Delta t)_{\text{meas}} - (\Delta t)_{\text{set}}$ vs time in fill



- Calorimeters from same laser have similar, but not necessarily identical $\Delta(\Delta t)$.

$(\Delta t)_{\text{meas}} - (\Delta t)_{\text{set}}$ vs time in fill




Conclusions

- The WFD timing can be evaluated both with InFill and STDP events.
- InFill pulses are more numerous, but less precise than STDP.
- Both show a stability within ~ 10 ps over the Fill (conclude from the errors).
- All 1296 crystals show a similar (and consistent) behaviour.

Back up Slides

Δt between consecutive pulses – Calo 1, all xtals

Δt between consecutive pulses (c.t)



t within fill pulses (μs)

Shows consistent behaviour of all xtals.

Bin size: Y – axis 0.1 ct, X- axis 2 μs
X, Y Ranges same as before

Δt = Difference between 2nd and 1st pulse for a fill

